

Virtuosity
is in your hands

X-SMART™

The endo motor for general dental practitioners



DENTSPLY
MAILLEFER

X-SMART™

compact light-weight simple use

Technical features

- Operates on battery or electrical power
- Battery recharging time: approx. 5 hours
- Time of use on battery: approx. 2 hours
- Supplied with a 16:1 contra-angle
- Operates without a pedal (sold separately)
- Also available combined to a *PROTAPER®* kit



Handpiece with integrated on/off switch (only 92g)

Large LCD screen

Autoreverse Mode

**Speed range*:
120 to 800 rpm**

**Torque control range*:
0,6 to 5,2 N.cm**

9 adjustable program selections

* with 16:1 contra-angle

X-SMART™

REF A 1004

X-SMART™ Contra-angle

REF A 1006

